











OTORHINOLARYNGOLOGY

- **High Quality** Accurate cutting edge made of high quality stainless steel. With improvement over previously available blades, these new blades combine quality and durability for sharp cutting.
- **Easy** These disposable blades are gamma sterilized and can be used directly after removal from the case.
- **Economical** After use, replace in the heat-resistant resin case for autoclaving or gas sterilization. The blade may also be re-used.



Myringotomy Blade

 <p>1.8mm</p>	 <p>75mm</p> <p>100% size</p> <p>K-18R</p>	Handle MF-100 MF-130 MF-130S
 <p>1.5mm</p>	 <p>69mm</p> <p>100% size</p> <p>K-15L</p> <p>K-15R</p>	
 <p>1.2mm</p>	 <p>63mm</p> <p>100% size</p> <p>K-12L</p> <p>K-12R</p>	

Laryngo Blade





 <p>1.5mm</p>	 <p>53.5mm</p> <p>100% size</p> <p>K-30</p>	Handle MF-200BA
--	---	---------------------------

Nasal Mucosa Blade

 <p>2.8mm</p> <p>10mm</p>	 <p>29.5mm</p> <p>100% size</p> <p>K-50S</p>	Handle MF-100 MF-130 MF-130S
--	--	--

Titanium Handle

The length shown below includes the base of the blade

MF-100	 <p>100mm*</p>
MF-130	 <p>130mm*</p>
MF-130S	 <p>130mm*</p>
MF-200BA	 <p>200mm*</p>



Mizuho America, Inc. P: 800.699.2547
 30057 Ahern Avenue F: 510.324.4545
 Union City, CA 94587 E: mizuho@mizuho.com
 mizuho.com



FEATHER SAFETY RAZOR CO., LTD.
 OVERSEAS TRADE DIVISION

MADE IN JAPAN