

REV. 120115



Distributed by:



Mizuho America, Inc.
30057 Ahern Avenue
Union City, CA 94587
(800) 699-CLIP (2547)
Fax: (510) 324-4545
mizuho.com

In partnership with:



REX



Single-Use, Non-Stick Bipolar Forceps



REX

SINGLE-USE, NON-STICK BIPOLAR FORCEPS

- Tips length 4.5 mm
- Shorter tips for precise coagulation
- Thinner and down angled bayonetted shaft for a better view under the microscope
- Ergonomic and knurled handle for a more secure grip
- Longer handle for better balance
- Compatible with 2-pin American Standard Connector
- Solid silver alloy with non-stick technology



CODE	L	T	QUANTITY
2606516	165 mm	0.25 mm	5
2606517	165 mm	0.5 mm	5
2606518	165 mm	1 mm	5
2606519	165 mm	1.5 mm	5
2606501	180 mm	0.25 mm	5
2606502	180 mm	0.5 mm	5
2606503	180 mm	1 mm	5
2606504	180 mm	1.5 mm	5
2606506	200 mm	0.25 mm	5
2606507	200 mm	0.5 mm	5
2606508	200 mm	1 mm	5
2606509	200 mm	1.5 mm	5
2606511	220 mm	0.25 mm	5
2606512	220 mm	0.5 mm	5
2606513	220 mm	1 mm	5
2606514	220 mm	1.5 mm	5

L = Forceps Length
T = Tip Width