

MST-7300 Mizuho Surgical Table Accessories


Standard Table Accessories



08-066-10-UE
Raised armboard (UE)
(08-117-09 required for use)




08-086-15-R1
Lower armboard (R1)



08-084-02
Body strap (1700 mm)



08-069-15-UE
Chest & waist support (UE)
(08-117-10 required for use)
(minimum 3 pcs required for use)




08-062-00-R1
Anesthesia screen
L type with clamp (R1)
** Wing type is also available.*




08-117-06 Standard armboard (C1)
08-117-07 60 mm Pad for armboard (C1/C2)
07-117-08 100 mm pad for armboard (C1/C2)
(minimum 2 pcs required for use)

Recommended Table Accessories



07-953-01
Sitting position frame
with side rail clamp
(08-117-09: 2 pcs required for use)



08-072-01
Knee rest
(08-117-09: 2 pcs required for use)



08-086-27
Arm shield with pad (130 mm)
*160 mm and 190 mm also available




08-066-05-NC
Armboard for prone position
(08-117-09 required for use)




08-069-24-NC
Lateral support
(08-117-10 required for use)
(minimum 2 pcs required for use)



08-080-11-NC
Adjustable body support
(08-117-09 required for use)
(minimum 2 pcs required for use)




08-070-25 Carbon fiber extension
08-070-32 70 mm Pad for Carbon
fiber extension



02-121-10
Radiolucent carbon fiber
four-point spinal frame



08-007-22-R1
Cervical Spine Frame (R1)



08-089-18
Radiographic top for 7300 series

MST-7300 Mizuho Surgical Table Accessories

Operation Table

Part No.	Description
18-350-70	Operating Table MST-7300B
18-351-70	Operating Table MST-7300BX

Standard Table Accessories

Part No.	Description	Quantity Required for Use	Additional Accessories	Compatibility (7201 VS 7300)
08-066-10-UE	Raised Armboard (UE)		08-117-09 (clamp for round bar)	Yes
08-086-15-R1	Lower Armboard (R1)			Yes
08-084-02	Body Strap (1700 mm)			Yes
08-069-15-UE	Chest and Waist Support (UE)	3	08-117-10 (clamp for flat bar)	Yes
08-062-00-R1	Anesthesia Screen L Type with Clamp (R1)			Yes
08-061-20-UE	Anesthesia Screen with Wing Type (UE)			Yes
08-117-06	Standard Armboard (C1)	2		Yes
08-117-07	60 mm Pad for Armboard (C1/C2)	2		Yes
08-117-08	100 mm Pad for Armboard (C1/C2)	2		Yes

Recommended Table Accessories

Part No.	Description	Quantity Required for Use	Additional Accessories	Compatibility (7201 VS 7300)
07-953-01	Sitting Position Frame With Side Rail Clamp		08-117-09 (clamp for round bar)	Yes
08-072-01	Knee Rest		08-117-09 (clamp for round bar)	Yes
08-086-27	Arm Shield with Pad (130 mm) <i>Includes two pcs</i>			Yes
08-066-05-NC	Armboard For Prone Position		08-117-09 (clamp for round bar)	Yes
08-069-24-NC	Lateral Support (NC)	2	08-117-10 (clamp for flat bar)	Yes
08-080-11-NC	Adjustable Body Support (NC)	2	08-117-09 (clamp for round bar)	Yes
08-070-25	Carbon Fiber Extension			
08-070-32	70 mm Pad For Carbon Fiber Extension			
02-121-10	Radiolucent Carbon Fiber Four-Point Spinal Frame			Yes
08-007-22-R1	Cervical Spine Frame (R1)			
08-089-18	Radiographic Top For 7300 Series			

TTA-7300-v01.1 Aug. 2017

Carefully read the Product Instructions prior to using products.

7201 and 7300 Series tables have different table-top dimensions:
Accessories for the 7300 Series table may not be compatible with the 7201 Series

Contact a Mizuho Sales Representative or Customer Service for more information.

CONTACT

Mizuho America, Inc.
30057 Ahern Avenue
Union City, CA 94587
www.mizuho.com
P: 800-699-CLIP (2547)

