

CERTIFICATE OF DECONTAMINATION

* This form must be filled in completely, signed, dated, and returned with the product(s) to Mizuho America, Inc. *

Hospital or Representative: _____

RGA#: _____

Description, Item Number, Serial and/or Lot Number:

OSHA regulation related to Bloodborne Pathogens (29 CFR, 1910.1030, including Subtitle B, Chapter XVII, Part 1910, Occupational Safety and Health Standards) prescribes safeguards to protect workers from contact with blood and other infectious waste. Additionally, the Department of Transportation regulations (49 CFR, 173.386 and 173.387) stipulates the safe transportation of etiologic agents.

In order to comply with Code of Federal Regulations and to ensure the safety of our employees and those that transport our product; Mizuho America requires that products sent for repair, returns, lab, or demo items used in a clinical environment, must be decontaminated. To provide objective evidence that these requirements have been met, this certificate must be completed prior to shipping.

PRODUCTS SHIPPED MUST NOT HAVE BEEN USED ON A PATIENT SUFFERING, OR SUSPECTED OF SUFFERING FROM VARIANT CREUTZFELDT-JAKOB DISEASE (vCJD).

Complete all applicable sections prior to shipping

- Has never been exposed in a clinical setting or environment
- Product is new and in its original packaging
- Product has been Decontaminated (reference AAMI Standard ST79)

A neutral pH detergent is preferred over caustic detergents. If using an enzyme cleaning immersion (neutral pH 6-8) detergent, a multi-step process is necessary to ensure thorough cleaning that may include thermal disinfection with an effective chelating detergent and mechanical actions to remove soils and microorganisms. Do not use alkaline, acid, or household detergents

Hospital Name, Phone Number & email: _____

Hospital Representative or MA Representative Name: _____

Title: _____ Signature: _____ Date: _____

The undersigned verifies that all products listed above have been properly processed according to the manufacturer's specifications prior to being shipped to Mizuho America, Inc.

Failure to comply with these policies and procedures may result in a minimum fine of \$1,500.00(USD).

For any questions or concerns contact Customer Service at 1-800-699-2547.